

SEQUENCHER[®]

Tutorial for Windows and Macintosh

Sequencher Connections

© 2015 Gene Codes Corporation

Gene Codes Corporation



Gene Codes Corporation
775 Technology Drive, Ann Arbor, MI 48108 USA
1.800.497.4939 (USA) +1.734.769.7249 (elsewhere)
+1.734.769.7074 (fax)
www.genecodes.com gcinfo@genecodes.com

Sequencher Connections

About Data Sources, Sessions, and Channels	3
Getting Started	3
Launching a Connections Session with BLAST	3
Viewing the Results	4
Changing Options for a BLAST Session	5
Launching a Connections Session with Primer-BLAST	6
Changing Options for your Primer-BLAST Session	7
Looking at the Results Using the Schematic	9
Launching a Connections Session with MUSCLE	10
Adding a Second Connections Session with MUSCLE	11
The Session is Saved with your Project	13
The Session Expires from NCBI	14
Restored Results	15
Deleting Sessions	15
Conclusion	15

Sequencher Connections

Sequencher Connections is a whole new way of performing multiple analyses on multiple sequences in parallel. This saves you time and gives you the opportunity to look at your data using different settings or parameters. Sequencher Connections, at the time of writing, allows you to use BLAST, Primer-BLAST, and Local-BLAST (if you have installed them) on individual sequences and MUSCLE on a group of related sequences. For these reasons we differentiate between individual and group or grouped data when talking about Connections. Not only can you view your BLAST, Primer-BLAST, Local-BLAST, or MUSCLE results at the time you performed the analysis, but you can also review your results later on since they are saved with your project. This tutorial covers working with sessions with BLAST, Primer-BLAST, and MUSCLE. For information on working with Local-BLAST in Connections, please refer to the Sequencher Connections chapter of the Sequencher User Manual or to the Local-BLAST tutorial.

ABOUT DATA SOURCES, SESSIONS, AND CHANNELS

A data source in the context of Sequencher Connections is a Sequencher project containing reads and contigs. A session is a coordinated set of data and analyses that operate on that data. A channel is an analysis such as MUSCLE or service such as BLAST. We will be using these terms throughout the tutorial.

GETTING STARTED

In this tutorial, you will learn how to use Sequencher Connections, how to create a session, how to name your session and change its channel options. You will also learn how to review a session that you created previously.


- Launch Sequencher.
- Go to the **File** menu and select **Import > Sequencher Project...**
- Navigate to the Sample Data folder inside the Sequencher application folder.
- Choose the **Sequencher Connections** project and select the **Open** button.
- Select **File > Save Project** and enter a name of **Sequencher Connections** and a location of your Desktop or Documents folder and select the **Save** button.
- Click on the **Label** column header in the Project Window to sort the project items by **Label**. This will make it easier to locate sequences later in the tutorial.

This project contains a set of fish mitochondrial sub-sequences. During the course of the tutorial, you will find out exactly what these sequences represent. These sequences will be analysed using BLAST. In the second part of the tutorial, you will analyse the sequences with Primer-BLAST and in part three, you will analyse the sequences with MUSCLE. There are also smaller sections on Restored Results and Deleting Sessions.


LAUNCHING A CONNECTIONS SESSION WITH BLAST

In this session you will use the default settings for BLAST. The sequences have been labeled using Sequencher's **Label** command to make it easier for you to choose the correct sequences in this tutorial.

- Select only the sequences labeled in blue (Tilapia). You should have selected 4 sequences.
- Go to the **Window** menu and select **Add to Connections Session...**
- Click on the radio button **Add to new Connections Session for individual sequences**.
- Type **Tilapia** in the **New Session Name** text box.



- Click the **OK** button.
- A new **Sequencher Connections** session starts. You will see a window with a table pane at the top and a pane at the bottom containing a series of tabs.
- Right-click on the **BLAST-nr** column header and choose **Run on Each Sequence**.
- Watch the table as the colors of each cell in the **BLAST-nr** channel column change from white to orange to cyan to yellow and then to green and the status finally changes to **Done**. Once all the cells contain a green **Done**, your analyses are complete and you may review the results.




- If you get any red **Failed** statuses in the cells, try to rerun the blasts for those particular cells by right-clicking in the cells and selecting **Run**.

VIEWING THE RESULTS

Once the BLAST run has finished, clicking in a **Done** cell will change the view in the Web View tab to the familiar BLAST results page. You can scroll up and down, and click on links.

- Click on the **Done** cell in the BLAST-nr column for the sequence called FISH067-08|Oreochromis in the Connections Session table.
- In the Web View tab, scroll down until you see the **Graphic Summary** representing the alignment of your query sequence and the alignments found by the BLAST algorithm.



- Click on the first red line. The web page jumps down to the alignment.

Download GenBank Graphics

Oreochromis niloticus voucher PNT-32 cytochrome oxidase subunit 1 (COI) gene, partial cds; mitochondrial
Sequence ID: [gb|KC789552.1](#) Length: 708 Number of Matches: 1

Range 1: 67 to 685

Score	Expect	Identities	Gaps	Strand
1144 bits(619)	0.0	619/619(100%)	0/619(0%)	Plus/Plus


Query 1	GGAAGTGCACCTAAGCCTCCTAATTCGGGCAGAACTAAGCCAGCCCGGCTCTCTCTCGGA	60
Sbjct 67	GGAAGTGCACCTAAGCCTCCTAATTCGGGCAGAACTAAGCCAGCCCGGCTCTCTCTCGGA	126
Query 61	GACGACCAATCTATAATGTAATTGTTACAGCACATGCTTTGTAATAATTTTCTTTATA	120
Sbjct 127	GACGACCAATCTATAATGTAATTGTTACAGCACATGCTTTGTAATAATTTTCTTTATA	186
Query 121	GTAATACCAATTATGATTGGAGGCTTTGGAACTGACTAGTACCCCTCATGATTGGTGCA	180
Sbjct 187	GTAATACCAATTATGATTGGAGGCTTTGGAACTGACTAGTACCCCTCATGATTGGTGCA	246

- Click on the **Sequence ID** link. A new window opens displaying the **Feature Table** for the hit sequence to which your target has aligned. Notice that this is mitochondrial Cytochrome oxidase subunit 1.
- Now click on a different green Done cell and see if you get the same or similar results for the other sequences.

CHANGING OPTIONS FOR A BLAST SESSION

BLAST offers a number of options for refining the scope of your search. Some of these options change values. One of the most important options is the E or expectation value. Another option you can change is the database to be searched. If you do not know much about the region you have just sequenced, you may want to cast the widest possible net by choosing the blastn algorithm rather than the megablast algorithm, which is the default setting. You can also ask BLAST to send you more than the default 10 results and descriptions as shown in this part of the tutorial below.

- Right-click in the header of the **BLAST-nr** channel.
- Choose **Insert BLAST Channel Before**.
- Right-click again in the header of the newly created **BLAST-nr** channel.
- Choose **Options...** from the context-sensitive menu.



- The **Channel Options** dialog opens. This dialog will contain a tab for each channel in your session.
- Click on **Optimize for** and change the algorithm to **Somewhat similar (blastn)**, clicking **OK** on the subsequent warning dialog.
- Change the word size to 11 to match the dialog you see below.
- Click on the **Default Graphic Color** picker, choose a new color for this channel (this will be visible in the Schematic), and click on the **OK** button.
- Change the number of alignments to 50 by choosing this value from the **Alignments** dropdown menu.
- Change the number of descriptions to 50 by choosing this value from the **Descriptions** dropdown menu.


- Click the **OK** button to save these changes. This closes the **Channel Options** dialog.
- Right-click on this channel's column header and choose **Run on Each Sequence**.
- If you get any red **Failed** statuses in the cells, try to rerun the blasts for those particular cells by right-clicking in the cells and selecting **Run**.

You will see more results and descriptions with the above settings than in the first part of the tutorial. Note that, although you see more results, these are still all mitochondrial fish sequences related to the initial sequences you chose as your queries. It will not always be the case that the results are so highly related, this will depend on your query sequence as well as which database you chose to search and the settings you used.

LAUNCHING A CONNECTIONS SESSION WITH PRIMER-BLAST

Primer-BLAST combines the capabilities of two resources into one. The first resource is the ability to harness the Primer 3 algorithm to predict primers and the second resource is the ability to use BLAST to check the specificity of those predictions. We will look at primer prediction only in this tutorial.

- Select the sequences labeled with magenta (Brittle star) in Sequencer's Project Window. You should have selected 4 sequences.
- Go to the **Window** menu and select **Add to Connections Session...**
- Click on the radio button **Add to new Connections Session for individual sequences**. Notice that your previous session is listed in the Session Launcher dialog.
- Type **Brittle star** into the **New Session Name** text box.



- Click the **OK** button.
- A new Sequencer Connections session starts. You will see a window with a table pane at the top and a pane at the bottom containing a series of tabs.
- Right-click on the **Primer-BLAST** column header and choose **Run on Each Sequence**.
- Watch the table as the colors of each cell change from white to cyan to green and the status changes to **Done**. Once all the cells contain a green **Done**, your analyses are complete and you may review the results.
- If you get any red **Failed** statuses in the cells, try to rerun the blasts for those particular cells by right-clicking in the cells and selecting **Run**.
- Click on the **Done** cell for the sequence called ACAP127-12|Ophiactis in the Primer-BLAST column. The Primer-BLAST results page loads in the **Web View** tab.

You can see your results displayed in a Graphical View with 10 sets of primer pairs in the display.

- Scroll down the page to view the detailed primer reports. You can see part of such a report in the image below.

Detailed primer reports

Primer pair 1								
	Sequence (5'→3')	Template strand	Length	Start	Stop	Tm	GC%	Self complementarity
Forward primer	AAGAAAGAGGAGCCGGAACG	Plus	20	138	157	60.04	55.00	4.00
Reverse primer	ACAGTGGCATTCCGTCCATT	Minus	20	335	316	59.96	50.00	3.00
Product length	198							

CHANGING OPTIONS FOR YOUR PRIMER-BLAST SESSION

Primer-BLAST offers a number of options for refining the scope of your search for primers. Some of these options are specific to the primer calculations but you may also check on the specificity of your primers. Note that if you do check your primers specificity, the search will take longer.

- Right-click in the header of the **Primer-BLAST** channel.

- Choose **Insert Primer-BLAST Channel Before**.
- Right-click in the header of the new **Primer-BLAST** channel and choose **Options...** from the context-sensitive menu.
- The **Channel Options** dialog opens. This dialog opens at the Primer-BLAST channel options tab for the new channel but will also contain a tab for each channel in your session.
- Click on the **Default Graphic Color** picker and choose a new color for this channel (this will be visible in the Schematic) and select the **OK** button.
- Change the **PCR product size** minimum from 70 to 500.
- Change the value of the **Number of primers to return** to 20.
- Change the **Max%** for **GC Content** to 75.

Primer-BLAST Search Tool Name: Primer-BLAST

This applies to **Individual Sequences** Default Graphic Color: [Blue Box]

Primer Parameters

Server location:


PCR product size: Minimum Maximum

Number of primers to return:

Primer melting temperatures: Minimum Optimum Maximum Maximum Tm Difference

GC Content: Min % Max % Clamp

- Click the **OK** button.
- Right-click on the new **Primer-BLAST** column header and choose **Run on Each Sequence**.
- Watch the table as the colors of each cell change from white to cyan to yellow to green and the status changes to **Done**. You should get a single red **Failed** status for one cell where Primer-BLAST failed to find any results within the minimum PCR product size which was the same as the length for this sequence.
- Click on one of the **Done** cells and click on the **Web View** tab.
- You can see your results displayed in a Graphical View with 20 sets of primer pairs in the display.




- Click on the **Text** tab. Now right-click in the tab view and choose **Save as Text...**, this saves a copy of the Primer-BLAST information to your chosen location.

Notice that if you insert the two sets of Primer-BLAST channels as neighbours, you can easily jump from one set of results to the other. In the next section, you will see how to view all the results for a single sequence using the Schematic.


LOOKING AT THE RESULTS USING THE SCHEMATIC

The Schematic enables you to combine several sets of different analyses into a single graphic. Thus for each sequence you are exploring, you can not only look at the results from several different BLAST searches but also combine that with the results of your Primer-BLAST analysis. This means that you can double check, with the aid of BLAST, that you are amplifying precisely the region you want and be assured of that in one single view instead of hopping about between different views or even different programs.

- Right-click on the sequence name in the second column of the session table and choose **Show Schematic**.



- A new dialog opens with the results of all the analyses you have run on your data in one view.
- Place your cursor over one of the colored lines in the **Schematic**, note the contents of the tooltip that appears. It contains information on the hit sequence in a **BLAST** search.
- Scroll in the window until you see **Primer-BLAST** results. The arrows represent the primers and the thin lines connecting the arrows represent the insert. Place your cursor over one of the colored lines and note the contents of the tooltip. It contains information on the primer pair.
- Click on any line in the schematic and see how it highlights a region corresponding to the range of the line you clicked on.
- Click on the radio button next to **Bases** and scroll to the bottom of the Schematic. You will see the bases contained within the range of the line you clicked on and their positions.



- In the bottom pane of the Schematic window, scroll down to view the table of results.

- Close the Schematic window.

LAUNCHING A CONNECTIONS SESSION WITH MUSCLE

In this session you will use the default settings for MUSCLE. The sequences have been labeled using Sequencer's **Label** command to make it easier for you to choose the correct sequences in this tutorial.

- Select the sequences labeled in cyan (Snapping shrimp). You should have selected 11 sequences.
- Go to the **Window** menu and select **Add to Connections Session...**
- Click on the radio button **Add to new Connections Session for grouped sequences**.
- Type **Snapping shrimp** into the **New Session Name** text box.

☐ Add to new Connections Session for individual sequences
☒ Add to new Connections Session for grouped sequences
 New Session Name:

- Click the **OK** button.

A new window opens containing the grouped sequences session. You will notice that there is only one line in the table and that it is called Group 1. The name is already highlighted so that you can type in a new, more memorable name. The table gives you information on how many sequences are contained in the group.

	Name	Num Seqs	MUSCLE-1	MUSCLE-2
1	Group 1	11	Queued	Queued


- Type the name **Synalpheus** to replace the words Group 1.

	Name	Num Seqs	MUSCLE-1	MUSCLE-2
1	Synalpheus	11	Queued	Queued

- Right-click on the MUSCLE-1 column header and choose **Run on Each Group**.
- The colors and text in the MUSCLE-1 column change in response to the status of the analysis.
- When the color in the cell is green and the text has changed to **Done**, click in the cell.
- The contents of the first tab change from the Connections landing page to a display of the aligned sequences.

Web View	Alignment	Log	Tree	Sequence
GBCD172-12_Synalpheus	-----CGGAATAATcTCTC			
GBCD168-12_Synalpheus	-----TCaC			
GBCD156-12_Synalpheus	-----TCCTAATtccCCnGCcITCGGtATAATcTCTC			
GBCD174-12_Synalpheus	-----TtCTAATcCTcCCgGCTTTCGGGAATAATcTCTC			
GBCD176-12_Synalpheus	-----TCCTAATTCTcCCaGCTTTCGGGAATAATTCcC			
GBCD155-12_Synalpheus	-----TtCTAATcTACCcGCcTtTGGGAATAATTCTC			
GBCD173-12_Synalpheus	-----TtTAATTCTnCCgGCTTTCGGATgATTCTC			
GBCD137-12_Synalpheus	gattctctgggacacctgaagtataTtTAATTCTnCCcGCTTTCGGATgATTCTC			
GBCD164-12_Synalpheus	-----CCTAATtTACCgGCTTTCGGtATAATTTCgC			
GBCD160-12_Synalpheus	-----tGGtATAATTTCaC			
GBCD165-12_Synalpheus	-----TtTAATTtTACCgGgTtTGGtATAATTTCaC			

- Click on the **Tree** tab.



- You now see a phylogram of the members of your data set.
- Right-click in the page and choose **Print Page...** Note that you could also change the tree from rectangular to a circular form or zoom in and out of the view (useful with large trees).
- Leave the session window open.

ADDING A SECOND CONNECTIONS SESSION WITH MUSCLE

When you add sequences to Sequencher Connections you can either add them to a new session or an existing session. If you are working with a set of related sequences you might want to add those to an existing group session so you can compare and contrast your results. The sequences are added as a new group and appear in the Sequencher Connections table as an additional row.

- Select all the sequences in the project - green (Sharks), magenta (Brittle stars), cyan (Snapping shrimp), and blue (Tilapia). You should have selected 24 sequences.
- Go to the **Window** menu and select **Add to Connections Session...**
- Click on the radio button **Add to an existing Connections Session.**
- Choose the **Snapping shrimp** group.

☒ Open an existing Connections window

Session Name	Type	Size	Created
Snapping shrimp	Individual	1	5/6/2015 3:19 PM
Brittle star	Individual	4	5/6/2015 3:19 PM
Tilapia	Individual	4	5/6/2015 3:05 PM

- Click the **OK** button.
- Notice that a new group has been added to the open session. Leave this name as **Group 2**.

	Name	Num Seqs	MUSCLE-1	MUSCLE-2
1	Synalpheus	11	Done	Queued
2	Group 2	24	Queued	Queued

- Right-click on the **MUSCLE-2** column header and choose **Run on Each Group**.

- When the color in both cells is green and the text has changed to **Done**, click in the Group 2 cell in the MUSCLE-2 channel. In the Web View, you will see the alignment for all the sequences.


Web View	Alignment	Log	Tree	Sequence
FISH005-08_Oreochromis	-----G			
FISH067-08_Oreochromis	-----G			
FISH052-08_Oreochromis	-----G			
TSM012-14_Oreochromis	-----atcGgcACCCTcTatCTAgTaTTTggTGCTTGAGCcGGAATAGTaG			
ANGBF1290-12_Carcharodon	-----CCTTTAttTAATTTTTggTGCaTGAGCAGGAATAGTgG			
FOAG009-07_Apristurus	-----CCTTaTACCTAATTTTTggTGCaTGAGCAGGcATgGtTg			
ENHHS111-06_Sphyrna	-----CCTTTACCTAATTTTTggTGCaTGAGCAGGtATAGTtG			
EWIGS001-06_Galeocerdo	-----CCTTTAtCTtATTTTTggTGCaTGAGCAGGtATAGTtG			
ESHKC083-07_Carcharhinus	-----CCTTTACCTgATTTTTggTGCaTGAGCAGGcATAGTtG			
EAR104-07_Ophiactis	-----cACCCTaTACtTtATTTTTggaGCTTGAGCAGGAacAGTaG			
ACAP127-12_Ophiactis	-----AACCCTTTAttTtAtaTTTggaGCaTGAGCTGGAacAGTaG			
EAR102-07_Ophiactis	-----CG			
GBMIN108-12_Ophiactis	-----			
GBCDA172-12_Synalpheus	-----			
GBCDA168-12_Synalpheus	-----tCCTTAATTTccCCnGCcT-----TCG			
GBCDA156-12_Synalpheus	-----tCCTAATTTcCCaGCTT-----TCG			
GBCDA176-12_Synalpheus	-----ttCTAATccTcCCgGCTT-----TCG			
GBCDA174-12_Synalpheus	-----tCtTAATcTtACcGCcT-----TtG			
GBCDA155-12_Synalpheus	-----tttTAATTCnCCgGCTT-----TCG			
GBCDA173-12_Synalpheus	gattctctgggacacctGaAggttaTattTAATTCnCCgGCTT-----TCG			
GBCDA137-12_Synalpheus	-----CCTAATTTtACaGCTT-----TCG			
GBCDA164-12_Synalpheus	-----tttTAATTTtACGgGgT-----TtG			
GBCDA165-12_Synalpheus	-----			
GBCDA160-12_Synalpheus	-----tG			

In this example, you can see that the Synalpheus portion of the alignment is much more fragmented, there are more gaps as a result of adding in the additional sequences, and fewer matching bases.

- Now switch to the **Alignment** tab. Here you scroll right and left to view the alignment rather than down and up with the results in the **Web View** tab.

Web View	Alignment	Log	Tree	Sequence
FISH005-08_Oreochromis	-----GGAAACCGCGCTAAGCCTCCTAATTCGGGCAGAA			
FISH067-08_Oreochromis	-----GGAACTGCACCTAAGCCTCCTAATTCGGGCAGAA			
FISH052-08_Oreochromis	-----GGAAACCGCACTAAGCCTCCTAATTCGGGCAGAA			
TSM012-14_Oreochromis	-----ATCGGCACCCCTCTATCTAGTATTGGTGCTTGAGCCGGAATAGTAGGAACCTGCAATTAAGCCTCCTAATTCGGGCAGAA			
ANGBF1290-12_Carcharodon	-----CCTTTATTTAATTTTTGGTGCTGAGCAGGAATAGTAGGAACAGCCCTTAAGCCTTTTAATCCGTGCCGAG			
FOAG009-07_Apristurus	-----CCTATACCTAATTTTTGGTGCTGAGCAGGATAGTTGGAATAGCCTTAAGTTTGTAAATCCGTGCCGAG			
ENHHS111-06_Sphyrna	-----CCTTTACCTAATTTTTGGTGCTGAGCAGGATAGTTGGAACAGCCCTTAAGCTTTTAATTCGAGCTGAA			
EWIGS001-06_Galeocerdo	-----CCTTTATCTTATTTTTGGTGCTGAGCAGGATAGTTGGAACAGCTCTAAGCTTTTAATTCGAGCTGAA			
ESHKC083-07_Carcharhinus	-----CCTTTACCTGATTTTTGGTGCTGAGCAGGATAGTTGGAACAGCTCTTACCTACTAATTCGAGCTGAA			
EAR104-07_Ophiactis	-----CACCCCTACTTTTATTTTTGGAGCTTGAGCAGGAACAGTAGGACTGCAATGAGAAAAATTATACGAGTCGAA			
ACAP127-12_Ophiactis	-----AACCCTTTATTTTATTTTTGGAGCATGAGCTGGAACAGTAGGCACATCCATGAGAAACATTATTTCGAGTAGAA			
EAR102-07_Ophiactis	-----AGTAGCC			
GBMIN108-12_Ophiactis	-----CGGAATAATCTCTCATATTATCAACCAAG-----TCACACATTATTAAACCAAG-----			
GBCDA172-12_Synalpheus	-----TCCTAATTTCCCCNGCCT-----TCGGTATAATCTCTCACATCATTAACCAAG-----			
GBCDA168-12_Synalpheus	-----TCCTAATTTCCCCAGCTT-----TCGGAATAATTTCCCAATCATTAACCAAG-----			
GBCDA156-12_Synalpheus	-----TTCTAATCTCTCCCGCTT-----TCGGAATAATCTCTCACATCATCAATCAAG-----			
GBCDA176-12_Synalpheus	-----TTCTAATCTTACCCGCTT-----TTGGAATAATTTCTCATATTATTAATCAAG-----			
GBCDA174-12_Synalpheus	-----TTTTAATCTCTNCCGCTT-----TCGGAATGATTTCTCATATTATTAATCAAG-----			
GBCDA155-12_Synalpheus	GATTCTTGGGACACCTGAAGTTTATTTTAATCTCTNCCGCTT-----TCGGAATGATTTCTCATATTATTAATCAAG-----			
GBCDA173-12_Synalpheus	-----CCTAATTTTACCAGCTT-----TCGGTATAATTCGCATATTATTAATCAAG-----			
GBCDA137-12_Synalpheus	-----TTTTAATTTTACCAGGCTT-----TTGGTATAATTCACATATTATTAATCAAG-----			
GBCDA164-12_Synalpheus	-----			
GBCDA165-12_Synalpheus	-----			
GBCDA160-12_Synalpheus	-----TGATATAATTCACATATTATTAATCAAG-----			

- Now switch to the **Tree** tab.



- Now click in the **Done** cell in the MUSCLE-2 Synalpheus group. In the **Tree** tab, you will see the tree for these sequences. In this way, you can visually compare the results of several sets of alignments quickly.

In this example, you added both sets of data to the same session, this makes sense for closely related data. If your data is less closely related, you can add it to a new Sequencher Connections session for grouped sequences instead.

THE SESSION IS SAVED WITH YOUR PROJECT

Your session is saved with your project and can be reviewed next time you open your project. You can test this with the following steps:

- Close the project.
- You should be prompted to save. Go ahead and Save.
- Exit Sequencher.

The next time you open the Sequencher Connections project you will be able to review the results from this tutorial as follows.

- Launch Sequencher.
- Go to the **File** menu and select **File > Open Project...**
- Navigate to the location you saved your project to earlier.
- Choose the **Sequencher Connections** project and select the **Open** button.
- Go to the **Window** menu and choose **Open Existing Connections Session...**
- **Save** the project if prompted to do so.

☐ Add to new Connections Session for individual sequences
☐ Add to new Connections Session for grouped sequences
 New Session Name:

☒ Open an existing Connections window

Session Name	Type	Size	Created
Snapping shrimp	Individual	1	5/6/2015 3:19 PM
Brittle star	Individual	4	5/6/2015 3:19 PM
Tilapia	Individual	4	5/6/2015 3:05 PM

- Choose the **Snapping shrimp** Grouped session by clicking on it. Notice that it shows 2 groups are associated with this session.
- Click the **OK** button.

The session opens and the saved results from the last time you used this project are now available for review.

THE SESSION EXPIRES FROM NCBI

If more than 36 hours have passed since you ran a BLAST or Primer-BLAST analysis, the results will have expired from the NCBI website. When you click on a Done status cell, Sequencer will display a different view in the Web View tab informing you that the results have expired. However, since Sequencer Connections has saved your results you will still be able to view them.

- Click on the **Text** tab to view the text version of your search results.
- Right-click on the sequence name in the second column of the session table and choose **Show Schematic** to view the graphic view of your results.

Web View
Text
XML
Sequence

Gene Codes Corporation

T C A G E N E
A G T C O D E S

Sequencer Connections

Your BLAST search results have expired from the NCBI website.

You can see a copy of your results in text format in the Text and Schematic views.

If you re-run your search, your currently saved results will be overwritten and you will see the new results from the BLAST website in this view.

If you want to run another search on your sequence, but also keep your old results, insert a channel to run your new search in. To insert the new channel, right-click in any channel header and choose either Insert BLAST Channel Before.

Note that any new search will also expire as indicated in the RID field of the returned results.

Sequencer Connections Individual: Tilapia

Name	Length	BLAST-rr	BLAST-rr	Primer-BLAST
1 FISH005-08...	619	Done	Done	Done
2 FISH052-08...	619	Done	Done	Done

Web View Text XML Sequence

BLAST

Basic Local Alignment Search Tool

Home Recent Results Saved Strategies Help

NCBI BLAST! blast suite! Formatting Results - 641TRRW13

Edit and Resubmit Save Search Strategies Formatting options Download

Nucleotide Sequence (619 letters)

Note that if you rerun any of the existing results and the runs complete, you will lose the originals. If some time has passed since you last ran a search and you suspect new data may have been added to the databases you searched against, add a new channel and just search using that channel, then use the Schematic view to compare the results.

RESTORED RESULTS

If your sessions have earlier search results associated with them, those results are restored when those sessions are reopened. Any previous green Done statuses will have gray Done statuses.

Likewise, if your sessions have earlier search results associated with them and you choose to run another query, then cancel that request before the query is complete, your earlier saved results are also restored.

DELETING SESSIONS

If you have accidentally added sessions or no longer need sessions, you can remove them from your project.

- Go to the **Window** menu and choose **Delete Existing Connections Session...**
- In the **Delete Session(s)** dialog, select the session or sessions you wish to delete and click on the **Delete Selected Session(s)** button. You will be prompted to **Cancel** or **Continue**.

If you select **Continue** and you have any open and/or running sessions, your running sessions will be stopped, closed, and removed from your project.

If you select **Continue** and you have no open and/or running sessions, your selected sessions will simply be removed from your project.

If you select **Cancel**, any running or open sessions will remain so and will not be removed from your project.

CONCLUSION

In this tutorial, you have worked with Sequencher Connections. You have learned how to add data to a session and how to give your session a meaningful name. You have also learned how to change session options. You have used BLAST, Primer-BLAST, and MUSCLE. You have learned how to view your results using the Schematic view and you have viewed your MUSCLE results as phylograms. You have also learned that your results are saved with your project and can be reviewed at a later time. Finally you have learned that although results expire from the NCBI website you can always view your results in the Text and Schematic views.